

FreeSurfer Multiprocessing Pipeline

Contents

Module freesurfer_wrapper	2
freesurfer_wrapper	2
Requirements	2
Preparation	2
Overview run	2
recon	2
edit	3
recon_edit	3
How to run recon	3
Create the input file	3
Run recon	3
How to run edit	4
Create the input file	4
Run edit	4
How to run recon_edit	4
Create the input file	4
Run recon_edit	4
How to check for completed runs and hard errors	5
Done	5
Error	5
How to restart after a computer failure	5
Quality control	5
Automated quality analysis	5
Manual quality analysis	5
Step 1	5
Step 2	6
Sub-modules	6
Module freesurfer_wrapper.run	6
Functions	6
Function argument_parser	6
Function edit	6
Function handle_workers	7
Function parse_input_file	7
Function recon	7
Function recon_edit	7
Function run_command	8
Function worker	8
Namespace freesurfer_wrapper.scripts	8
Sub-modules	8
Module freesurfer_wrapper.scripts.check_logs	8
Functions	8
Function get_logs	8
Function print_id_from_logs	9

Module freesurfer_wrapper.scripts.create_recon_input	9
Functions	9
Function create_input_file	9
Module freesurfer_wrapper.scripts.view	9
Functions	10
Function freeview	10
Function get_subjects	10

Module freesurfer_wrapper

freesurfer_wrapper

freesurfer_wrapper aims to facilitate the creation of a multiprocessing pipeline using FreeSurfer. It is a Python wrapper to execute parallel runs of recon-all and some pial edits algorithms.

Requirements

- Docker¹
- FreeSurfer license key²

Preparation

- 1) Place the license in a license.txt file in the same folder as the Dockerfile.
- 2) Place your dataset folder in the same folder as the Dockerfile.
- 3) Build the docker image:

```
docker build -t fs_wrapper .
```

Overview run

The main script has 3 commands, each of them is explained below.

```
python run.py -h
usage: run.py [-h] {recon,edit,recon_edit} ...
```

Command-line wrapper tool to execute parallel runs of FreeSurfer recon-all and some pial edits algorithms.

positional arguments:

```
{recon,edit,recon_edit}
    recon          Run FreeSurfer recon-all.
    edit           Run mri_gcut and mri_binarize for pial edits.
    recon_edit     Re-run recon-all for pial edits.
```

optional arguments:

```
-h, --help          show this help message and exit
```

recon

Run FreeSurfer recon-all.

```
python run.py recon -h
usage: run.py recon [-h] -i INPUT [-p PARALLEL]
```

optional arguments:

¹<https://www.docker.com/>

²<https://surfer.nmr.mgh.harvard.edu/registration.html>

```

-h, --help            show this help message and exit
-i INPUT, --input INPUT
                        Tab separated file. First column: unique ID. Second
                        column: path to dcm/nii file.
-p PARALLEL, --parallel PARALLEL
                        Number of parallel runs (default: number of CPUs).

```

edit

Run mri_gcut and mri_binarize for pial edits.

```

python run.py edit -h
usage: run.py edit [-h] -i INPUT [-p PARALLEL]

```

optional arguments:

```

-h, --help            show this help message and exit
-i INPUT, --input INPUT
                        Tab separated file. First column: unique ID. Second
                        column: path to dcm/nii file. Third column: tissue
                        ratio
-p PARALLEL, --parallel PARALLEL
                        Number of parallel FS runs (default: number of CPUs).

```

recon_edit

Re-run recon-all for pial edits.

```

python run.py recon_edit -h
usage: run.py recon_edit [-h] -i INPUT [-p PARALLEL]

```

optional arguments:

```

-h, --help            show this help message and exit
-i INPUT, --input INPUT
                        Subject id list file.
-p PARALLEL, --parallel PARALLEL
                        Number of parallel FS runs (default: number of CPUs).

```

How to run recon

Create the input file

A tab-separated file is needed as input with the following characteristics:

- Each line must represent a single scan.
- First column: unique ID.
- Second column: path to dcm/nii file.

For the ADNI dataset, you can create this file using the create_recon_input.py script. This will create a recon_input.txt file. The script will combine the subject ID and the session ID to create a unique ID.

```

docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python scripts/create_recon_input.py

```

As an example, a recon_input.txt file for some ADNI records will look like this:

```

137_S_1414_S46193    ADNI/137_S_1414/MP-RAGE/2008-02-26_11_57_53.0/S46193/ADNI_137_S_1414_MR_MP-RAGE_
137_S_1414_S72806    ADNI/137_S_1414/MP-RAGE/2009-08-26_11_06_33.0/S72806/ADNI_137_S_1414_MR_MP-RAGE_

```

For other datasets you can try to edit the PATH_PATTERN variable in scripts/create_recon_input.py.

Run recon

Now you can run the recon command using Docker. This can take several hours.

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python run.py recon -i recon_input.txt
```

This command will use the maximum number of CPUs. You can append the `-p <INT>` flag where `<INT>` is the number of parallel runs you want.

How to run edit

Create the input file

After running recon you can check your results using freeview. Please refer to [Manual quality analysis section](#) to use a custom script.

If any skull edits are necessary, you need to create an input table to run edit with the following characteristics:

- Each line must represent a single scan.
- First column: unique ID.
- Second column: path to dcm/nii file.
- Third column: the tissue ratio for WM edits.

Tissue ratio is the threshold to value (%) of WM intensity. The value should be >0 and <1 ; larger values would correspond to cleaner skull-strip but higher chance of brain erosion.

You can copy the recon input table, keep only the lines for the scans that need pial edits, and add the column with the tissue ratio values.

As an example, a `edit_input.txt` file for some ADNI records will look like this:

```
137_S_1414_S46193    ADNI/137_S_1414/MP-RAGE/2008-02-26_11_57_53.0/S46193/ADNI_137_S_1414_MR_MP-RAGE_
137_S_1414_S72806    ADNI/137_S_1414/MP-RAGE/2009-08-26_11_06_33.0/S72806/ADNI_137_S_1414_MR_MP-RAGE_
```

Run edit

Now you can run the edit command using Docker.

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python run.py edit -i edit_input.txt
```

This command will use the maximum number of CPUs. You can append the `-p <INT>` flag where `<INT>` is the number of parallel runs you want.

You can check the resulting edited mask using freeview:

```
SUBJECTS_DIR=$(pwd)/FS_OUTPUTS
freeview -recon <UNIQUE_ID> -v brainmask.gcutsT$<TISSUE_RATIO>.mgz:colormap=heat:opacity=0.5
```

If still not good, change the tissue ratio value in the input file and run edit again. When all masks are OK, proceed to recon_edit command.

How to run recon_edit

Create the input file

The input file is the table used for edit with the final values for tissue ratio.

Run recon_edit

recon_edit will re-run parts of FS recon-all using the edited masks. This can take several hours.

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper python run.py recon_edit -i edit_i
```

This command will use the maximum number of CPUs. You can append the `-p <INT>` flag where `<INT>` is the number of parallel runs you want.

How to check for completed runs and hard errors

FreeSurfer's recon-all command creates different logs while running. The recon-all.done log is created only for completed runs. The recon-all.error is created for hard failures.

You can check these logs using a custom script. The script was written to work on ADNI folder structure. For other datasets you can try to edit the PATH_PATTERN variable in scripts/check_logs.py.

Done

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python scripts/check_logs.py done
```

You can also pipe the output to bash word count command to get a quick count:

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python scripts/check_logs.py done | wc -l
```

Error

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python scripts/check_logs.py error
```

How to restart after a computer failure

If the execution of the pipeline is halted by a computer failure or system restart then you have to update the input file of the recon command.

```
bash scripts/update_recon_input.sh
```

This will remove all “done” samples from the original recon_input. It will also delete the folders from the samples that were running when the failure happened. These samples will run again from scratch.

Now run again the [recon command](#) but using the new input file (<YYYY-MM-DD>_recon_input.txt).

Quality control

Automated quality analysis

The tool is packaged with **qatools-python**³ version 1.2 for quality control measurements. This script was developed by Reuter DeepMI Lab⁴ as a revision, extension, and translation to the Python language of the FreeSurfer QA Tools.

```
docker run --rm -it -v $(pwd):/root/freesurfer_wrapper fs_wrapper \
python scripts/qatools-python/qatools.py --subjects_dir FS_OUTPUTS --output_dir QC \
--screenshots --outlier --fornix
```

This will create qatools-results.csv file; screenshots, outliers and fornix folders inside the QC folder. Please consult qatools-python docs⁵ for a full explanation of each QC measurement.

Manual quality analysis

You can visually inspect each result using freeview. We provide a script to speed up the opening process of each scan. The script also prompts the user about the result of the QC after each window of freeview is closed. The result is saved to manual_QC.txt

It is not possible to run freeview using Docker, the graphical user interface cannot be displayed. Therefore, you need to have FreeSurfer/freeview installed in your host machine.

Step 1

Set SUBJECTS_DIR environment variable. Here, the results are stored inside the FS_OUTPUTS directory.

```
export SUBJECTS_DIR=$(pwd)/FS_OUTPUTS
```

³<https://github.com/Deep-MI/qatools-python>

⁴<https://deep-mi.org/>

⁵[scripts/qatools-python/README.md#description](#)

Step 2

Run the view.py script.

```
python scripts/view.py
```

Sub-modules

- [freesurfer_wrapper.run](#)
- [freesurfer_wrapper.scripts](#)

Module freesurfer_wrapper.run

Command-line wrapper tool to execute parallel runs of FreeSurfer recon-all and some pial edits algorithms.

This file can also be imported as a module and contains the following functions:

```
* argument_parser - parser for command-line options, arguments and sub-commands.
* run_command -
* handle_workers - creates a pool of parallel worker processes running commands.
* worker - invokes a subprocess running the command.
* recon - formats recon-all command string.
* edit - formats mri_gcut and mri_binarize command string.
* recon_edit - formats a cp and recon-all command string.
* parse_input_file - parses the input tables.
```

Functions

Function argument_parser

```
def argument_parser(
    args: list
) -> ArgumentParser.parse_args
```

Parser for command-line options, arguments and sub-commands.

Parameters

args : list Command-line arguments list

Returns

Parser

Function edit

```
def edit(
    edit_args: list
) -> str
```

Formats mri_gcut and mri_binarize command string. mri_gcut performs skull stripping algorithm based on graph cut. mri_binarize binarizes the edited mask.

Parameters

edit_args : list mri_gcut and mri_binarize arguments list

Returns

```
mri_gcut [args] && mri_binarize [args]
```

Function `handle_workers`

```
def handle_workers(  
    p: int,  
    command: function,  
    input_file: str  
)
```

Creates a pool of parallel worker processes running commands. Workers will be called until all lines from the input file are processed.

Parameters

p : int The number of parallel processes.
command : function Function returning the command-line string to pass the worker.
input_file : str Tab-separated .txt file.

Returns

None

Function `parse_input_file`

```
def parse_input_file(  
    input_file: str  
) -> List[List[str]]
```

Parses the input tables.

Parameters

input_file : str Tab-separated .txt file.

Returns

File lines and columns parsed as a list of lists.

Function `recon`

```
def recon(  
    recon_args: list  
) -> str
```

Formats recon-all command string.

Parameters

recon_args : list recon-all arguments list

Returns

recon-all [args]

Function `recon_edit`

```
def recon_edit(  
    recon_edit_args: list  
) -> str
```

Formats a cp and recon-all command string. cp replaces the original brainmask with the edited brainmask.gcutsT{tissue_ratio}.mgz. recon-all re-runs -autorecon2-wm -autorecon3 stream with the new mask.

Parameters

recon_edit_args : list cp and recon-all arguments list

Returns

cp [args] && recon-all [args]

Function `run_command`

```
def run_command(  
    args  
)
```

Pass the appropriate command function to the worker handler.

Parameters

args : list Command-line arguments list

Returns

None

Function `worker`

```
def worker(  
    cmd: str  
) -> <function run at 0x7f3aa69e6670>
```

Invokes a subprocess running the command.

Parameters

cmd : str Command-line string

Returns

`subprocess.run()`

Namespace `freesurfer_wrapper.scripts`

Sub-modules

- [freesurfer_wrapper.scripts.check_logs](#)
- [freesurfer_wrapper.scripts.create_recon_input](#)
- [freesurfer_wrapper.scripts.view](#)

Module `freesurfer_wrapper.scripts.check_logs`

Script to check recon-all logs for each run

usage: python check_logs.py <done|error>

Please edit the **PATH_PATTERN** variable with the appropriate pathname pattern to find each file.

This file can also be imported as a module and contains the following functions:

- * `get_logs` - get the path for each log based on pathname pattern.
- * `print_id_from_logs` - prints the IDs from a list of logs.

Functions

Function `get_logs`

```
def get_logs(  
    path_pattern: str  
) -> list
```

Get the path for each log based on pathname pattern.

Parameters

path_pattern : str Glob pathname pattern to find each log.

Returns

List of log paths

Function print_id_from_logs

```
def print_id_from_logs(  
    logs: list  
)
```

Prints the IDs from a list of logs.

Parameters

logs : list List of log paths

Returns

None

Module freesurfer_wrapper.scripts.create_recon_input

Script to create recon input table

This script creates an input table based on the directory organization of the image files.

Please edit the **PATH_PATTERN** variable with the appropriate pathname pattern to find each file.

This file can also be imported as a module and contains the following functions:

* create_input_file - creates the input table.

Functions

Function create_input_file

```
def create_input_file(  
    path_pattern: str  
)
```

Creates a two column text file to be used as input for the main script recon command. First column: unique ID (combines SUBJECT ID and SESSION ID). Second column: path to DICOM file.

Parameters

path_pattern : str Glob pathname pattern to find each DICOM file.

Returns

None

Module freesurfer_wrapper.scripts.view

Script to open the scans in sequence and register if passed or not in QC analysis

This script identifies the subjects present in the folder set in SUBJECTS_DIR environment variable. Using this list, it automatically opens freeview and asks the user input for the result of QC analysis. Results are saved in manual_QC.txt file.

This file can also be imported as a module and contains the following functions:

* get_subjects - returns a list of subjects inside the folder.
* freeview - returns a freeview command formatted as string

Functions

Function freeview

```
def freeview(  
    subject_id: str  
)
```

Returns a freeview command formatted as string.

Parameters

subject_id : str FreeSurfer subject_id.

Returns

str

Function get_subjects

```
def get_subjects(  
    subjects_dir='/subjects'  
)
```

Returns a list of subjects inside the folder.

Parameters

subjects_dir : str, default=os.environ['SUBJECTS_DIR'] FreeSurfer SUBJECTS_DIR.

Returns

list

Generated by *pdoc* 0.10.0 (<https://pdoc3.github.io>).